

**Prescription Ware Ordering Guide**

HSAH Item #	Quantity per Case	Centor Item #	Size
-------------	-------------------	---------------	------

	<b>Bulk Pack- Screw-Loc® Amber Plastic Vials &amp; Closures</b> (Child resistant closure imprinted with HSAH logo)			
	018169	475 / cs	D-6-BTR	6 dram
	018170	355 / cs	D-8-BTR	8½ dram
	018171	250 / cs	D-13-BTR	13 dram
	018172	200 / cs	D-16-BTR	16 dram
	018173	150 / cs	D-20-BTR	20 dram
	018174	120 / cs	D-30-BTR	30 dram
	018175	85 / cs	D-40-BTR	40 dram
	018176	55 / cs	D-60-BTR	60 dram

	<b>Bulk Pack- Screw-Loc® Green Plastic Vials &amp; Closures</b> (Child resistant closure imprinted with HSAH logo)			
	024211	475 / cs	DG-6-BTR	6 dram
	024212	355 / cs	DG-8-BTR	8½ dram
	024213	250 / cs	DG-13-BTR	13 dram
	024214	200 / cs	DG-16-BTR	16 dram
	024215	150 / cs	DG-20-BTR	20 dram
	024216	120 / cs	DG-30-BTR	30 dram
	024217	85 / cs	DG-40-BTR	40 dram
	024218	55 / cs	DG-60-BTR	60 dram


	<b>Bulk Pack- Screw-Loc® Blue Plastic Vials &amp; Closures</b> (Child resistant closure imprinted with generic logo)			
	038958	475 / cs	DX-6-BSH	6 dram
	038959	355 / cs	DX-8-BSH	8½ dram
	038960	250 / cs	DX-13-BSH	13 dram
	038961	200 / cs	DX-16-BSH	16 dram
	038962	150 / cs	DX-20-BSH	20 dram
	038964	120 / cs	DX-30-BSH	30 dram
	038965	85 / cs	DX-40-BSH	40 dram
	038966	55 / cs	DX-60-BSH	60 dram

	<b>Amber Glass Dropper Service with Plastic Pipettes</b> (Black molded cap and rubber bulb)			
	004356	12 / carton	0764-07A	½ oz.
	004357	12 / carton	0765-07A	1 oz.
	004358	12 / carton	0766-06A	2 oz.


HSAH Item #	Quantity per Case	Centor Item #	Size
-------------	-------------------	---------------	------

	<b>Bulk Pack- 1-Clic® Amber Plastic Vials &amp; Closures</b> (Reversible, white closure- imprinted with HSAH logo)			
	018162	200 / cs	B-10-BTR	10 dram
	018163	175 / cs	B-13-BTR	13 dram
	018164	150 / cs	B-16-BTR	16 dram
	018165	120 / cs	B-20-BTR	20 dram
	018166	80 / cs	B-30-BTR	30 dram
	018167	70 / cs	B-40-BTR	40 dram
	018168	45 / cs	B-60-BTR	60 dram

	<b>Bulk Pack- 1-Clic® Green Plastic Vials &amp; Closures</b> (Reversible, white closure- imprinted with HSAH logo)			
	024219	200 / cs	BG-10-BTR	10 dram
	024220	175 / cs	BG-13-BTR	13 dram
	024221	150 / cs	BG-16-BTR	16 dram
	024222	120 / cs	BG-20-BTR	20 dram
	024223	80 / cs	BG-30-BTR	30 dram
	024224	70 / cs	BG-40-BTR	40 dram
	024225	45 / cs	BG-60-BTR	60 dram

	<b>Amber Plastic Ovals - Graduated</b> (White child resistant Clic-Loc III® caps attached)			
	032918	132 / cs	PE-2A	2 oz.
	032919	110 / cs	PE-3A	3 oz.
	032920	90 / cs	PE-4A	4 oz.
	032921	72 / cs	PE-6A	6 oz.
	004874	50 / cs	PE-8	8 oz.
	032922	42 / cs	PE-12A	12 oz.
032923	36 / cs	PE-16A	16 oz.	

	<b>Amber Glass Ovals - Graduated</b> (White child resistant Clic-Loc III® caps attached)			
	004249	72 / cs	P-801	1 oz.
	004250	72 / cs	P-804	4 oz.

	<b>White Plastic Ointment Jars</b>			
	032925	12 / bag	PL-3701A	1 oz.
	032926	12 / bag	PL-3704A	4 oz.
	032927	16 / bag	PL-3708A	8 oz.
	032928	9 / bag	PL-3716A	16 oz.

