



# KatKath<sup>TM</sup>

TomCat Urinary Catheter  
*Be kind to your patients*

## Benefits

- Easy placement
- Flexible for maximum comfort
- Thin walled-larger bore lumen
- Great longitudinal strength
- Rounded atraumatic tip



Rotating Suture Collar



Increased Flexibility

### More Flexible/Low coefficient of friction

The catheter is made from a unique non-reactive material which has similar properties to Telfon® but is more flexible. More flexibility allows greater comfort.

### Crush Resistant

The KatKath<sup>TM</sup> is almost impossible to crush and is also highly resistant to kinking. This allows placement without the use of a stylet. (Useful if flushing is required during the act of placement).

### Well Tolerated

The KatKath<sup>TM</sup> may be left in place for hours or days, dependant on individual needs, and is extremely well tolerated by even the most difficult patients.

### Great Longitudinal Strength

The KatKath's increased flexibility does not diminish its ability in the effective treatment of blocked cats due to its great longitudinal strength and will not concertina even when encountering gross resistance.

### Rotating Suture Collar

The KatKath<sup>TM</sup> has a rotating suture collar-To allow natural rotation while in situ, this also lowers tensions on sutures.


Order Code	Details	Open Ended
J-1203	3.5Fr OD. 11cm long with Luer lock	
J-1203A	3.5Fr OD. 14cm long with Luer lock	

**JorVet™**


# KatKath™

**Adjustable Urinary Catheter**

*Be kind to your patients*



Adjustable hub slides to desired length



Suture groove

## Benefits

**Easy placement**  
**Flexible for maximum comfort**  
**Thin walled-larger bore lumen**  
**Great longitudinal strength**  
**Rounded atraumatic tip**

## Features

**Non-reactive material**  
**Suture groove on hub**  
**Low coefficient of friction**  
**Kink resistant**  
**Crush resistant**  
**Lightweight**  
**Radiopaque**  
**No stylet required**

### **Order Code**

J-1203B

### **Details    Open Ended**

3.5Fr OD. Up to 18cm maximum usable length with Luer lock

Jorgensen Laboratories

1450 Van Buren Ave

Loveland, CO 80538

1-800-525-5614

www.JorVet.com    Info@JorVet.com