



KatKathTM

TomCat Urinary Catheter
Be kind to your patients

Benefits

- Easy placement
- Flexible for maximum comfort
- Thin walled-larger bore lumen
- Great longitudinal strength
- Rounded atraumatic tip



Rotating Suture Collar



Increased Flexibility

More Flexible/Low coefficient of friction

The catheter is made from a unique non-reactive material which has similar properties to Telfon® but is more flexible. More flexibility allows greater comfort.

Crush Resistant

The KatKathTM is almost impossible to crush and is also highly resistant to kinking. This allows placement without the use of a stylet. (Useful if flushing is required during the act of placement).

Well Tolerated

The KatKathTM may be left in place for hours or days, dependant on individual needs, and is extremely well tolerated by even the most difficult patients.

Great Longitudinal Strength

The KatKath's increased flexibility does not diminish its ability in the effective treatment of blocked cats due to its great longitudinal strength and will not concertina even when encountering gross resistance.

Rotating Suture Collar

The KatKathTM has a rotating suture collar-To allow natural rotation while in situ, this also lowers tensions on sutures.


Order Code	Details	Open Ended
J-1203	3.5Fr OD. 11cm long with Luer lock	
J-1203A	3.5Fr OD. 14cm long with Luer lock	

JorVet™


KatKath™

Adjustable Urinary Catheter

Be kind to your patients



Adjustable hub slides to desired length



Suture groove

Benefits

Easy placement
Flexible for maximum comfort
Thin walled-larger bore lumen
Great longitudinal strength
Rounded atraumatic tip

Features

Non-reactive material
Suture groove on hub
Low coefficient of friction
Kink resistant
Crush resistant
Lightweight
Radiopaque
No stylet required

Order Code
J-1203B

Details Open Ended
3.5Fr OD. Up to 18cm maximum
usable length with Luer lock

Jorgensen Laboratories
1450 Van Buren Ave
Loveland, CO 80538
1-800-525-5614
www.JorVet.com Info@JorVet.com