



## Medical Indications:

Diets with these nutritional modifications have been recommended for cats with the following conditions:

- Diabetes mellitus
- Persistent hyperglycemia
- Enteritis, diarrhea
- Hepatic lipidosis
- Critical care management

## Medical Contraindications:

- Renal failure
- Hepatic encephalopathy



Actual Size



**AVAILABLE SIZES**  
Bags of 6 and 10 lbs  
Cases of 24 5.5-oz cans

## Feeding and Administration:

Animal feeding tests using Association of American Feed Control Officials (AAFCO) procedures substantiate that Purina Pro Plan Veterinary Diets DM Dietetic Management Feline Formulas provide complete and balanced nutrition for maintenance of adult cats.

The following feeding program is recommended as a guideline only, with discretionary clinical adjustments for proper weight maintenance.

FEEDING GUIDE <sup>1</sup>			
BODY WEIGHT		DRY	CANNED
lbs	kg	Approximate Number Of 8-oz Cups/Day <sup>3</sup>	Approximate Number Of 5.5-oz Cans/Day <sup>3</sup>
5	2.3	1/4	3/4
8	3.6	1/3	1 1/8
12	5.4	1/2	1 3/4
16	7.2	3/4	2 1/4
20	9.1	1	2 7/8
25 <sup>4</sup>	11.4	1 1/8	3 1/2
CALORIE CONTENT (FED) <sup>2</sup>			
		4118 kcal/kg	1227 kcal/kg
		1868 kcal/lb	557 kcal/lb
		592 kcal/8-oz cup	191 kcal/5.5-oz can

## Glucose Management:



High protein and low carbohydrate to help support the unique nutritional needs of diabetic cats



May help promote insulin sensitivity and reduce external insulin requirements



May help reduce post-prandial hyperglycemia

## General Health:

- St/Ox Defense (dry and wet formulas): This diet promotes a urinary environment unfavorable to the development of both struvite and calcium oxalate crystals
- High protein levels help maintain lean body mass and ideal body condition
- High in antioxidant Vitamins A and E to support a healthy immune system
- Also available in Savory Selects to help with long-term compliance
- Great taste to help ensure adequate nutrient intake

## Ingredients (Dry):

Poultry by-product meal, soy protein isolate, corn gluten meal, soy flour, animal fat preserved with mixed-tocopherols, corn starch, animal liver flavor, calcium carbonate, phosphoric acid, fish oil, potassium chloride, L-Lysine monohydrochloride, DL-Methionine, taurine, choline chloride, powdered cellulose, salt, zinc sulfate, Vitamin E supplement, ferrous sulfate, manganese sulfate, niacin, Vitamin A supplement, calcium pantothenate, copper sulfate, thiamine mononitrate, riboflavin supplement, Vitamin B-12 supplement, pyridoxine hydrochloride, folic acid, Vitamin D-3 supplement, calcium iodate, biotin, menadione sodium bisulfite complex (source of Vitamin K activity), sodium selenite. A-2561



## Ingredients (Canned):

Liver, poultry by-products, meat by-products, water sufficient for processing, chicken, salmon, oat fiber, salmon meal, artificial and natural flavors, calcium sulfate, guar gum, potassium chloride, carrageenan, salt, Vitamin E supplement, mono and dicalcium phosphate, taurine, thiamine mononitrate, zinc sulfate, ferrous sulfate, niacin, calcium pantothenate, copper sulfate, Vitamin A supplement, manganese sulfate, menadione sodium bisulfite complex (source of Vitamin K activity), pyridoxine hydrochloride, riboflavin supplement, Vitamin B-12 supplement, biotin, folic acid, Vitamin D-3 supplement, potassium iodide. A-2530



### AVERAGE NUTRIENT CONTENT

	DRY			CANNED		
	PER 100 KCAL ME	AS FED	DRY MATTER	PER 100 KCAL ME	AS FED	DRY MATTER
Protein	12.75 g	52.50%	56.97%	10.76 g	13.20%	53.93%
Fat	4.02 g	16.55%	17.96%	5.34 g	6.55%	26.76%
Carbohydrate	3.58 g	14.75%	16.01%	1.44 g	1.77%	7.23%
Crude Fiber	0.26 g	1.09%	1.18%	0.60 g	0.74%	3.02%
Total Dietary Fiber	2.67 g	11.00%	11.94%	2.45 g	3.01%	12.30%
Soluble Fiber	0.10 g	0.40%	0.43%	0.23 g	0.28%	1.14%
Insoluble Fiber	2.57 g	10.60%	11.50%	2.22 g	2.73%	11.15%
Calcium	0.26 g	1.09%	1.18%	0.27 g	0.33%	1.35%
Phosphorus	0.26 g	1.09%	1.18%	0.26 g	0.32%	1.31%
Potassium	0.21 g	0.85%	0.92%	0.20 g	0.25%	1.02%
Sodium	0.12 g	0.51%	0.55%	0.10 g	0.12%	0.49%
Chloride	0.17 g	0.69%	0.74%	0.15 g	0.19%	0.78%
Magnesium	0.02 g	0.10%	0.10%	0.02 g	0.02%	0.08%
Copper	0.50 mg	20.6 mg/kg	22.4 mg/kg	0.39 mg	4.8 mg/kg	19.6 mg/kg
Zinc	5.95 mg	245 mg/kg	263 mg/kg	4.32 mg	53 mg/kg	217 mg/kg
Taurine	0.08 g	0.32%	0.35%	0.06 g	0.07%	0.29
Vitamin A	695 IU	28,620 IU/kg	31,059 IU/kg	4,980 IU	61,101 IU/kg	249,653 IU/kg
Vitamin E	12.94 IU	533 IU/kg	578 IU/kg	7.91 IU	97 IU/kg	396 IU/kg
Total Omega-6 Fatty Acids	0.43 g	1.78%	1.93%	0.86 g	1.06%	4.33%
Total Omega-3 Fatty Acids	0.08 g	0.32%	0.35%	0.10 g	0.12%	0.48%
EPA+DHA	0.05 g	0.19%	0.20%	0.06 g	0.07%	0.30%

### DIGESTION TEST RESULTS†

DIGESTIBILITY:	DRY	CANNED
Total, %	85.3	81.4
Protein, %	92.3	88.5
Fat, %	90.1	89.2
Carbohydrate, %	84.4	30.3
Calorie, %	89.4	85.1
PERCENTAGE OF METABOLIZABLE ENERGY FROM:		
Protein, %	48.9	42.8
Fat, %	37.4	51.5
Carbohydrate, %	13.7	5.7

†Based on digestion testing conducted at the Purina PetCare Technology Centers

1 Adjustments must be made for environmental conditions, activity level, and body condition. Provide fresh water in a clean container daily.

2 Kilocalories of metabolizable energy (ME)

3 This daily amount should be divided into several small meals throughout the day.

4 For each additional pound of body weight, feed an additional 24 kilocalories.