



## JorVet Warming Air Blankets

Keeping patients warm during surgery and recovery has always been a challenge. Electric pads have a risk of overheating, circulating water blankets are easily damaged.

Our new warm air blankets are a safe, economical solution to this problem. They are made of a silky polyester fabric with an inlay of nylon water proof material to keep the airflow circulating upward towards the patient.

These blankets are machine washable (hang to dry) therefore much more economical than other air blankets on the market, which are disposable. There is a 5" diameter round opening with a nylon cord to secure the blanket to any of the major brands of warm air convection units. For Veterinary use only.

Blanket sizes available:

<b>J1034</b>	Small	30"L x 20"W
<b>J1034A</b>	Medium	36"L x 24"W
<b>J1034B</b>	Large	48"L x 24"W
<b>J1034S</b>	All 3	

Air is filtered via 0.2 micron filter inside the convection unit.



Warm air flows only upward

### Disposable Blanket

One side has nylon reflective coating, the other (patient) side is SMS material. Light blue. 27" x 34"

### J1034F

Will attach to all brands of warmers.

Also available through Jorgensen Laboratories the **J-1034C** Warm Air Unit by CSZ.

Features:

- 3 temperature settings
- Audio / visual alarms
- 6 ft. flexible hose
- Dimensions: 8.75"W x 8.75"D x 13.5"H
- 13.5 lbs.
- 110 - 120V, 50-60Hz
- Filtration - 0.2 microns, High Efficiency

### Jorgensen Laboratories, Inc.

1450 North Van Buren Avenue  
Loveland, CO 80538  
(800) 525-5614  
www.JorVet.com