

NEW

AniSet IV Fluid Sets



Features Include:

- Compatible with all IV accessories from the major brands
- Designed and configured to work with Abbott, Heska and Hospira IV pumps that are calibrated with Abbott settings.
- Compare & Save! Extremely competitive pricing

Anti-kink tube: allows free movement of your patient without kinking or blockage.

Set Break Point: includes extension set and allows patients to be moved without the need to remove or re-prime the set.

Auto-Fill Air Cap: allows the set to automatically prime and fill without spilling out the end. No more messes or wasted fluid!

With
**AniSet the
Extension
set is
FREE!**



Aniset Anti-Kink IV 78" L/L (15 drps/ml)
Aniset Anti-Kink IV 78" S/L (15 drps/ml)
Aniset Anti-Kink IV 106" L/L (15 drps/ml)
Aniset Anti-Kink IV Pediatric/Microbore 70" L/L (60 drps/ml)
Aniset IV Opset 73" S/L (15 drps/ml)
Male Adapter Plug S/L (pk 1)

**Anti-Kink
Tubing**

The only IV set designed with animals in mind

BROUGHT TO YOU BY
CLIPPER
www.clipperdist.net