

Cell Cycle 360

Bench Top Centrifuge

Medivet Cell Cycle 360 offers precise digital timing and speed control, combined with easy to navigate touch panel and highly visible LED display.


Ideal for veterinary laboratory centrifugation, the Cell Cycle 360 has the option of up to 4 interchangeable rotor/ carrier combinations suitable for tubes ranging from 5ml to 50ml.

- Brushless DC motor
- Digital LED display
- Soft touch control panel
- Electronic lid interlock with alarm
- Stainless steel bowl
- Interchangeable swing and fixed angle rotors
- Adjustable acceleration and deceleration
- RCF and rpm display
- Veterinary Use Only



Cell Cycle 360

TECHNICAL SPECIFICATION	
Catalogue	CEN 91720
Maximum Speed *	5,000 rpm
Maximum RCF *	4,100 g
Timer	0 - 99 min
Power Rating	120 V 60 Hz 5A
End of Cycle	Audible alert
Dimensions	250 mm H 432 mm D 370 mm W
Weight	15 kg

ROTORS					
Catalogue Number	Description	Maximum speed [rpm]	Maximum RCF [g]	Capacity	
CEN 91720	Cell Cycle 360 Centrifuge (120V 60Hz) (rotor not included)				
CEN 91700 004	Swing rotor for Cell Cycle 360 complete with 4 x 50ml stainless carriers	4,200	2,760	4 x 50 ml	

ACCESSORIES	
Catalogue Number	Description
CEN 91720 017	User Manual for CEN 91720 Cell Cycle 360 Centrifuge
CEN 91710 034	Adaptor, 50ml - 15ml Falcon™ / 8.5 ml Vacutainer® Tube (Pack 4)
CEN 91700 035	Adaptor, 50ml - 6ml Vacuette® Tube (Pack 4)
CEN 91700 036	Tool, Rotor Screwdriver, Phillips 110 x 5.5mm
CEN 91700 037	Tool, Lid Release Screwdriver, Phillips 120 x 2.9mm
CEN 91700 039	Lead, Mains Power, IEC to Earthed American Plug

All specifications are subject to change without notice. Please check the exact specification when ordering.
 * Subject to rotor selection

Medivet America INC
 200 Moore Drive
 Nicholasville KY 40356

T: 1-800-365-9168
 E: info@medivet-america.com
 W: www.medivet-america.com

