



COOLRENEWAL
www.cool-renewal.com

*Providing a Low Cost,
High Value
Cryosurgical Solution
to Veterinarians
around the World*



Benefits

- Minimal Discomfort
- Self-Anesthetizing
- 2 Minute Procedure
- No Suturing
- No Evaporation
- No Expiration
- Portable
- Low Cost Treatment
- 3 Disposable Applicator Options
- Multiple Applicator Sizes for Various Lesion Sizes and Locations

Commonly Treated Skin Lesions

Papillomas



Sebaceous Adenoma



Skin Tag




Distichia



Histiocytoma



Sarcoid



Other treatable lesions include, but are not limited to:

Follicular Cysts, Granulomas, Fibroids, & more!

Cool Renewal™ is a convenient and cost effective alternative to traditional liquid nitrogen methods and equipment, as well as other expensive portable cryosurgery systems. Cool Renewal™ was developed to fill the gaps of its competitor products, offering veterinarians a wide variety of applicators in several sizes, new applicators unique to our cryosurgical system, and the lowest cost per treatment available from a portable cryosurgical system.

Cool Renewal™ Complete Cryosurgical Kit

Item # J1401



- ❖ (2) 170mL Canisters of Cryogen
~ 130 Freezes*
 - ❖ (2) Extender Tubes
 - ❖ Assorted Pack of 60 Disposable, Double Ended Foam Tipped Applicators
(20 ea, Small, Medium, Large/XL)
 - ❖ Assorted Pack of 50 Disposable Plastic Isolation Funnels
(10 ea- 5mm, 8mm, 12mm, 14mm, 16mm)
 - ❖ (20) Disposable Skin Tag Tweezers
 - ❖ Instructions for Use
 - ❖ Patient Instructions Tear Pad
 - ❖ Applicator Practice Pad
 - ❖ Magnetic Packaging for Component Storage
- *Average Number of Freezes based on treatment of 5mm lesions

Cool Renewal™ Replacement Canister

Item # J1401B



- ❖ Single Replacement Canister-
170ml (65 Freezes)
- ❖ Extender Tube
- ❖ Instructions for Use
- ❖ Patient Instructions Tear Pad

Foam Tipped Applicators



- ❖ Small- **Item# J1401D1**
- ❖ Medium- **Item# J1401D2**
- ❖ Large/ XL- **Item# J1401D3**
- ❖ Assorted- **Item# J1401D4**
*60 per pack

Skin Tag Tweezers



- ❖ **Item # J1401A**
(qty 20)

Isolation Funnels



- ❖ 5mm- **Item#J1401D5**
(qty 50)
- ❖ 8mm- **Item#J1401D6**
(qty 50)
- ❖ 12mm- **Item#J1401D7**
(qty 50)
- ❖ 14mm- **Item#J1401D8**
(qty 50)
- ❖ 16mm- **Item#J1401D9**
(qty 50)
- ❖ Assorted-
Item#J1401D10
(10 of each size)

All applicators are disposable. DO NOT ATTEMPT TO CLEAN OR USE ON MULTIPLE PATIENTS