

# Start Saving Today with the Pureline® OC8200 Oxygen Concentrator.

**Pureline**®  
OXYGEN CONCENTRATORS

**Stop paying for tanks, refills, and deliveries. Stop wasting time on tank billing, scheduling, management and maintenance. Instead, get your oxygen for free from room air!**



## Savings for you, safety for your patients.

Supera has perfected the technology by engineering in safety and durability features not found on traditional oxygen concentrators.

All Pureline concentrators are designed specifically for use with veterinary anesthesia machines and oxygen therapy.

The built-in E-tank manifold allows mobile use of the machine. It also protects your patient in case of a power failure.

### Key Features

- Delivers USP93% oxygen up to 8 LPM
- Connects with industry standard oxygen DISS fittings
- Safety backup oxygen supply
- Energy efficient; Just pennies a day to operate
- Built specifically for anesthesia machine use
- Oxygen purity and “no flow” safety alarms
- Made in the USA by craftsmen using state of the art manufacturing equipment
- Virtually maintenance free

**The M8200 pays for itself in as little as one year.**

**Visit [www.superavet.com/ROI](http://www.superavet.com/ROI) to see how!**



# The Pureline® OC8200 Pays For Itself Fast!

*By eliminating the cost of oxygen tank rental, purchase, refills, delivery charges and tank management time, the Pureline OC8200 will pay for itself very quickly.*



## Return On Investment

can be in just months!



Visit [www.superavet.com/ROI](http://www.superavet.com/ROI) to calculate  
YOUR actual savings!



The heavy duty shelf is pre-drilled with VESA 100 pattern for easy mounting of monitors or other equipment.



Dual industry standard demand DISS connection. The auxiliary and main supply switch is clearly marked and accessible.

## The Benefits of the Pureline OC8200

- No tanks, hoses, regulators to buy, store, or replace
- Oxygen is "on demand"—never run out!
- Quiet 53-dBA operation
- Stable and strong "A" frame chassis with a low center of gravity
- Small 18" square footprint saves valuable floor space
- Includes one high pressure "E" tank for backup and oxygen flush (the back up e-tank is not filled by machine)

**SUPERA** | ANESTHESIA INNOVATIONS

Formerly LEI Medical

**Anesthesia Innovations Since 1993**  
(503) 723-5068 or (877) 620-1500 toll-free  
info@superavet.com  
12330 S.E. Highway 212, Clackamas, OR 97015