

EN GASTROENTERIC

LOW FAT™

CANINE FORMULA

Medical Indications:

Diets with these nutritional modifications have been recommended for dogs with the following conditions:

- Hyperlipidemia
- Lymphangiectasia
- Pancreatitis
- Protein-losing enteropathy
- Inflammatory bowel disease
- Exocrine pancreatic insufficiency
- Vomiting and diarrhea

Medical Contraindications:

- None



Actual Size

AVAILABLE SIZES

Bags of 6, 18, and 32 lbs
Cases of 12 13.4-oz cans



Feeding and Administration:

Animal feeding tests using Association of American Feed Control Officials (AAFCO) procedures substantiate that Purina Pro Plan Veterinary Diets EN Gastroenteric Low Fat Dry Canine Formula provides complete and balanced nutrition for maintenance of adult dogs. Purina Pro Plan Veterinary Diets EN Gastroenteric Low Fat Canned Canine Formula is formulated to meet the nutritional levels established by AAFCO for maintenance of adult dogs.

The following feeding program is recommended as a guideline only, with discretionary clinical adjustments for proper weight maintenance.

FEEDING GUIDE ¹			
BODY WEIGHT		DRY	CANNED
lbs	kg	Approximate Number Of 8-oz Cups/Day ³	Approximate Number Of 13.4-oz Cans/Day ³
5	2.3	3/4	2/3
10	4.5	1 1/8	1
20	9.1	1 7/8	1 2/3
30	13.6	2 1/2	2 1/4
40	18.2	3	2 2/3
50	22.7	3 1/2	3 1/8
60	27.3	3 7/8	3 1/4
80	36.4	4 3/4	4 1/4
100 ⁴	45.5	5 1/2	4 7/8
CALORIE CONTENT (FED) ²			
		3550 kcal/kg	959 kcal/kg
		1610 kcal/lb	435 kcal/lb
		327 kcal/8-oz cup	364 kcal/13.4-oz can

Gastrointestinal Health:



High total digestibility to promote nutrient absorption



Low in fat for dogs who have difficulty digesting fat



Contains the prebiotic inulin, soluble fiber that helps nourish beneficial microflora to promote a healthy GI tract

General Health:


- Added antioxidant Vitamins A and E to support a healthy immune system, including a high level of Vitamin E help to reduce oxidative stress
- Great taste to help ensure adequate nutrient intake

Ingredients (Dry):

Brewers rice, barley, corn gluten meal, poultry by-product meal, animal digest, animal fat preserved with mixed-tocopherols, calcium carbonate, mono and dicalcium phosphate, inulin, potassium chloride, sodium bicarbonate, salt, L-Lysine monohydrochloride, Vitamin E supplement, zinc proteinate, choline chloride, L-ascorbyl-2-polyphosphate (source of Vitamin C), manganese proteinate, ferrous sulfate, niacin, copper proteinate, Vitamin A supplement, calcium pantothenate, thiamine mononitrate, Vitamin B-12 supplement, riboflavin supplement, pyridoxine hydrochloride, garlic oil, folic acid, calcium iodate, menadione sodium bisulfite complex (source of Vitamin K activity), Vitamin D-3 supplement, biotin, sodium selenite. B-2633

Ingredients (Canned):

Water sufficient for processing, rice, liver, meat by-products, barley, soy protein isolate, potato protein, turkey, salmon, glycerin, powdered cellulose, guar gum, inulin, calcium sulfate, tricalcium phosphate, potassium chloride, carrageenan, magnesium sulfate, Vitamin E supplement, calcium carbonate, L-Lysine monohydrochloride, zinc sulfate, niacin, ferrous sulfate, thiamine mononitrate, natural smoke flavor, calcium pantothenate, pyridoxine hydrochloride, copper sulfate, Vitamin B-12 supplement, riboflavin supplement, Vitamin A supplement, manganese sulfate, folic acid, Vitamin D-3 supplement, biotin, potassium iodide, sodium selenite. B-2589



AVERAGE NUTRIENT CONTENT

	DRY			CANNED		
	Per 100 kcal ME	As Fed	Dry Matter	Per 100 kcal ME	As Fed	Dry Matter
Protein	7.32 g	26.00%	28.44%	12.17 g	11.67%	44.04%
Fat	1.74 g	6.19%	6.77%	2.38 g	2.28%	8.60%
Carbohydrate	14.28 g	50.70%	55.46%	10.31 g	9.89%	37.32%
Crude Fiber	0.31 g	1.10%	1.20%	0.96 g	0.92%	3.47%
Total Dietary Fiber	1.79 g	6.36%	6.96%	2.96 g	2.84%	10.72%
Soluble Fiber	0.32 g	1.13%	1.24%	0.67 g	0.64%	2.42%
Insoluble Fiber	1.47 g	5.23%	5.72%	2.29 g	2.20%	8.30%
Calcium	0.33 g	1.17%	1.28%	0.32 g	0.31%	1.17%
Phosphorus	0.23 g	0.80%	0.88%	0.27 g	0.26%	0.98%
Potassium	0.21 g	0.73%	0.80%	0.24 g	0.23%	0.87%
Sodium	0.14 g	0.48%	0.53%	0.15 g	0.14%	0.53%
Chloride	0.25 g	0.90%	0.98%	0.20 g	0.19%	0.72%
Magnesium	0.03 g	0.09%	0.10%	0.03 g	0.02%	0.09%
Copper	0.41 mg	14.6 mg/kg	16.0 mg/kg	0.43 mg	4.1 mg/kg	15.5 mg/kg
Zinc	6.20 mg	220 mg/kg	241 mg/kg	3.75 mg	36 mg/kg	136 mg/kg
Vitamin A	654 IU	23,206 IU/kg	25,386 IU/kg	2,985 IU	28,622 IU/kg	108,008 IU/kg
Vitamin E	18.37 IU	652 IU/kg	713 IU/kg	13.76 IU	132 IU/kg	498 IU/kg
Total Omega-6 Fatty Acids	0.42 g	1.48%	1.62%	0.46 g	0.45%	1.68%
Total Omega-3 Fatty Acids	0.02 g	0.07%	0.07%	0.05 g	0.04%	0.17%
EPA+DHA	0.00 g	0.00%	0.00%	0.02 g	0.02%	0.08%

DIGESTION TEST RESULTS†

DIGESTIBILITY:	DRY	CANNED
Total, %	87.8	81.1
Protein, %	87.1	79.5
Fat, %	86.7	84.0
Carbohydrate, %	94.5	91.2
Calorie, %	89.5	81.9
PERCENTAGE OF METABOLIZABLE ENERGY FROM:		
Protein, %	28.3	43.1
Fat, %	16.4	20.4
Carbohydrate, %	55.3	36.5

†Based on digestion testing conducted at the Purina PetCare Technology Centers

1 Adjustments must be made for environmental conditions, activity level, body condition and breed size. As reflected in the feeding tables, smaller breed dogs can require 15% to 25% more kilocalories per pound of body weight than larger breed dogs. Provide fresh water in a clean container daily.

2 Kilocalories of metabolizable energy (ME)

3 This daily amount should be divided into several small meals throughout the day.

4 For each additional 10 lbs of body weight, feed an additional 114 kilocalories.