

Designed specifically for veterinary applications and available in a broad range of sizes, Henry Schein Brand hard pack disposable needles, syringes and combinations consistently perform at high levels for smooth tissue penetration and improved patient comfort.

HARD PACK



Recyclable, hard "rigid" packaging is autoclavable, provides protection against contamination and additional safeguard during product transport.



Needles

- Sharp tri-bevel point and inert lubrication for smooth tissue penetration
- ISO recognized color-coding for accurate, quick size identification
- Aluminum insert for maximum stability
- Sterilized by ethylene oxide
- Universal design works with all Luer Lock or Luer Slip syringes

Syringes:

- Available with Luer Lock and Luer Slip tips
- Catheter and eccentric tips available on select sizes
- Polypropylene barrel and elastomer piston for smooth plunger motion
- Bold barrel markings for easy and accurate dosage measurements
- Barrel designed to prevent accidental plunger removal
- Large flanges for improved control
- Tipped plunger head on 1mL size minimizes dead space
- Sterilized by ethylene oxide
- Piston is free of natural latex

KEY FEATURES

Tri-Bevel Needle



Aluminum Insert



Bold Barrel Markings



Large Flanges



Quality, Reliability and Savings

Needles



SKU#	BARREL	NEEDLE GAUGE	LUER	UOM
062379	—	● 18 X 1	—	100 Count
062380	—	● 18 X 1-1/2	—	100 Count
062381	—	● 19 X 1	—	100 Count
062382	—	● 19 X 1-1/2	—	100 Count
062384	—	● 20 X 1	—	100 Count
062385	—	● 20 X 1-1/2	—	100 Count
062383	—	● 20 X 3/4	—	100 Count
062386	—	● 21 X 1	—	100 Count
062388	—	● 22 X 1	—	100 Count
062389	—	● 22 X 1-1/2	—	100 Count
062387	—	● 22 X 3/4	—	100 Count
062391	—	● 23 X 1	—	100 Count
062390	—	● 23 X 3/4	—	100 Count
062392	—	● 25 X 5/8	—	100 Count

Syringes



SKU#	BARREL	NEEDLE GAUGE	LUER	UOM
062396	3mL	—	Luer Lock	100 Count
062409	6mL	—	Luer Lock	50 Count
062411	12mL	—	Luer Lock	80 Count
062413	20mL	—	Luer Lock	50 Count
062416	35mL	—	Luer Lock	30 Count
062419	60mL	—	Luer Lock	20 Count
062394	1mL	—	Luer Slip	100 Count
062397	3mL	—	Luer Slip	100 Count
062410	6mL	—	Luer Slip	50 Count
062412	12mL	—	Luer Slip	80 Count
062414	20mL	—	Luer Slip	50 Count
062417	35mL	—	Luer Slip	30 Count
062420	60mL	—	Luer Slip	20 Count
062418	60mL	—	Catheter Tip	20 Count

Needle/Syringe Combinations



SKU#	BARREL	NEEDLE GAUGE	LUER	UOM
062398	3mL	● 20 X 3/4	Luer Lock	100 Count
062399	3mL	● 20 X 1	Luer Lock	100 Count
062401	3mL	● 20 X 1-1/2	Luer Lock	100 Count
062403	3mL	● 21 X 1	Luer Lock	100 Count
062404	3mL	● 22 X 3/4	Luer Lock	100 Count
062406	3mL	● 22 X 1	Luer Lock	100 Count
062408	3mL	● 25 X 5/8	Luer Lock	100 Count
062400	3mL	● 20 X 1	Luer Slip	100 Count
062402	3mL	● 20 X 1-1/2	Luer Slip	100 Count
062405	3mL	● 22 X 3/4	Luer Slip	100 Count
062407	3mL	● 22 X 1	Luer Slip	100 Count
062395	1mL	● 25 X 5/8	Luer Slip	100 Count

© 2018 Henry Schein Animal Health. All Rights Reserved. No copying without permission. Not responsible for typographical errors. Lit. No. 19HSB5 (01/2019)



Quality, Reliability and Savings