

VENDOR	SKU	DESCRIPTION	UOM	PACK TYPE	ACTIVE INGREDIENTS	WITHDRAWAL
DRY						
Boehringer Ingelheim	034456	Tomorrow Infusion 10mL Tube	144 Count	Pail	Cephapirin Benzathine (300mg)	42 days-meat/ 72 hrs post calving
Boehringer Ingelheim	034081	Tomorrow Mastitis Treatment 10mL Syringe	12 Count	Box	Cephapirin Benzathine (300mg)	42 days-meat/ 72 hrs post calving
Boehringer Ingelheim	004122	Dry Clox Infusion	144 Syringe	Pail	Cloxacillin Benzathine (500mg)	30 days-meat/ 30 days pre calving
Merck	021923	Orbenin-DC Dry Cow Intramammary Infusion	12 x 7.5gm	Box	Benazthine Cloxacillin (500mg)	28 days-meat/ 28 days pre calving
WG Critical Care	050474	Quartermaster Suspension	144 x 10mL	Pail	Procaine Penicillin G (1,000,000 units), Dihydrostreptomycin (1gm)	60 days- meat/ 96 hrs milk
WG Critical Care	050475	Quartermaster Suspension	12 x 10mL	Box	Procaine Penicillin G (1,000,000 units), Dihydrostreptomycin (1gm)	60 days- meat/ 96 hrs milk
Zoetis	019993	Albadry Plus Mastitis Tubes	144 x 10mL	Box	Penicillin G Procaine (200,000 IU) Novibocin Sodium (400mg), Chlorobutanol Anhydrous (50mg)	30 days-meat/ 72 hrs post calving
Zoetis	034200	Orbeseal Mastitis Tubes	144 Count	Pail	65% Bismuth Subnitrate (4gm)	0 days- meat/ 0 days milk
Zoetis	034117	Spectramast DC	12 Count	Syringe	Ceftiofur Equivalents (500mg), Microcrystalline Wax (700mg), Labrafil M 1944CS (500mg)	16 days-meat/ 0 days milk
LACTATING						
Boehringer Ingelheim	034455	Today Infusion 10ml Tube	144 Count	Pail	Cephapirin Sodium (200mg)	4 days -meat/ 96 hours milk
Boehringer Ingelheim	034080	Today Mastitis Treatment 10mL Syringe	12 Count	Box	Cephapirin Sodium (200mg)	4 days -meat/ 96 hours milk
Boehringer Ingelheim	058503	Polymast Infusion 10mL Syringe	12 Count	Box	Hetacillin Potassium	10 days-meat/ 72 hrs milk
Boehringer Ingelheim	058508	Polymast Infusion 10mL Syringe	144 Count	Box	Hetacillin Potassium	10 days-meat/ 72 hrs milk
Merck	021925	Amoximast Intramammary Infusion	12 x 10mL Syringe	Box	Amoxicillin (62.5mg)	12 days meat/ 60 hrs milk
Merck	021924	Dariclox Intramammary Infusion	12 x 10mL Syringe	Box	Sodium Cloxacillin (200mg)	10 days meat/ 48 hrs milk
Zoetis	034113	Pirsue Sterile Solution	12 Count	Syringe	Pirlimycin Hydrochloride (50mg)	9 days-meat/ 36 hrs milk
Zoetis	034114	Pirsue Sterile Solution Syringe	144 Count	Box	Pirlimycin Hydrochloride (50mg)	9 days-meat/ 36 hrs milk
Zoetis	034115	Spectramast LC	12 Count	Syringe	Ceftiofur Equivalents (125mg), Microcrystalline Wax (700mg), Labrafil M 1944CS (500mg)	2 days meat/ 72 hrs milk
Zoetis	034116	Spectramast LC Syringe	144 Count	Box	Ceftiofur Equivalents (125mg), Microcrystalline Wax (700mg), Labrafil M 1944CS (500mg)	2 days meat/ 72 hrs milk
Zoetis	034118	Spectramast DC Syringe	144 Count	Box	Ceftiofur Equivalents (500mg), Microcrystalline Wax (700mg), Labrafil M 1944CS (500mg)	16 days-meat/ 0 days milk



Mastitis tube reference chart



VENDOR	SKU	DESCRIPTION	UOM	PACK TYPE	ACTIVE INGREDIENTS	WITHDRAWAL
DRY						
Boehringer Ingelheim	034456	Tomorrow Infusion 10mL Tube	144 Count	Pail	Cephapirin Benzathine (300MG)	42 days-meat/ 72 hrs post calving
Boehringer Ingelheim	034081	Tomorrow Mastitis Treatment 10mL Syringe	12 Count	Box	Cephapirin Benzathine (300MG)	42 days-meat/ 72 hrs post calving
Boehringer Ingelheim	004122	Dry Clox Infusion	144 SYR	Pail	Cloxacillin Benzathine (500MG)	30 days-meat/ 30 days pre calving
Merck	021923	Orbenin-DC Dry Cow Intramammary Infusion	12 x 7.5gm	Box	Benazthine Cloxacillin (500MG)	28 days-meat/ 28 days pre calving
WG Critical Care	050474	Quartermaster Suspension	144 x 10mL	Pail	Procaine Penicillin G (1,000,000 units), Dihydrostreptomycin (1GM)	60 days- meat/ 96 hrs milk
WG Critical Care	050475	Quartermaster Suspension	12 x 10mL	Box	Procaine Penicillin G (1,000,000 units), Dihydrostreptomycin (1GM)	60 days- meat/ 96 hrs milk
Zoetis	019993	Albadry Plus Mastitis Tubes	144 x 10mL	Box	Penicillin G Procaine (200,000 IU) Novibiocin Sodium (400MG), Chlorobutanol Anhydrous (50MG)	30 days-meat/ 72 hrs post calving
Zoetis	034200	Orbeseal Mastitis Tubes	144 Count	Pail	65% Bismuth Subnitrate (4G)	0 days- meat/ 0 days milk
Zoetis	034117	Spectramast DC	12 Count	Syringe	Ceftiofur Equivalents (500MG), Microcrystalline Wax (700MG), Labrafil M 1944CS (500MG)	16 days-meat/ 0 days milk
LACTATING						
Boehringer Ingelheim	034455	Today Infusion 10ml Tube	144 Count	Pail	Cephapirin Sodium (200MG)	4 days -meat/ 96 hours milk
Boehringer Ingelheim	034080	Today Mastitis Treatment 10mL Syringe	12 Count	Box	Cephapirin Sodium (200MG)	4 days -meat/ 96 hours milk
Boehringer Ingelheim	058503	Polymast Infusion 10mL Syringe	12 Count	Box	Hetacillin Potassium	10 days-meat/ 72 hrs milk
Boehringer Ingelheim	058508	Polymast Infusion 10mL Syringe	144 Count	Box	Hetacillin Potassium	10 days-meat/ 72 hrs milk
Merck	021925	Amoximast Intramammary Infusion	12 x 10mL Syringe	Box	Amoxicillin (62.5MG)	12 days meat/ 60 hrs milk
Merck	021924	Dariclox Intramammary Infusion	12 x 10mL Syringe	Box	Sodium Cloxacillin (200MG)	10 days meat/ 48 hrs milk
Zoetis	034113	Pirsue Sterile Solution	12 Count	Syringe	Pirlimycin Hydrochloride (50MG)	9 days-meat/ 36 hrs milk
Zoetis	034114	Pirsue Sterile Solution Syringe	144 Count	Box	Pirlimycin Hydrochloride (50MG)	9 days-meat/ 36 hrs milk
Zoetis	034115	Spectramast LC	12 Count	Syringe	Ceftiofur Equivalents (125MG), Microcrystalline Wax (700MG), Labrafil M 1944CS (500MG)	2 days meat/ 72 hrs milk
Zoetis	034116	Spectramast LC Syringe	144 Count	Box	Ceftiofur Equivalents (125MG), Microcrystalline Wax (700MG), Labrafil M 1944CS (500MG)	2 days meat/ 72 hrs milk
Zoetis	034118	Spectramast DC Syringe	144 Count	Box	Ceftiofur Equivalents (500MG), Microcrystalline Wax (700MG), Labrafil M 1944CS (500MG)	16 days-meat/ 0 days milk

VENDOR	SKU	DESCRIPTION	UOM	PACK TYPE	ACTIVE INGREDIENTS	DRY/LACTATING	WITHDRAWAL
Boehringer Ingelheim	034455	Today Infusion 10ml Tube	144 Count	Pail	Cephapirin Sodium (200MG)	Lactating	4 days -meat/ 96 hours milk
Boehringer Ingelheim	034080	Today Mastitis Treatment 10mL Syringe	12 Count	Box	Cephapirin Sodium (200MG)	Lactating	4 days -meat/ 96 hours milk
Boehringer Ingelheim	058503	Polymast Infusion 10mL Syringe	12 Count	Box	Hetacillin Potassium	Lactating	10 days-meat/ 72 hrs milk
Boehringer Ingelheim	058508	Polymast Infusion 10mL Syringe	144 Count	Box	Hetacillin Potassium	Lactating	10 days-meat/ 72 hrs milk
Merck	021925	Amoximast Intramammary Infusion	12 x 10mL Syringe	Box	Amoxicillin (62.5MG)	Lactating	12 days meat/ 60 hrs milk
Merck	021924	Dariclox Intramammary Infusion	12 x 10mL Syringe	Box	Sodium Cloxacillin (200MG)	Lactating	10 days meat/ 48 hrs milk
Zoetis	034113	Pirsue Sterile Solution	12 Count	Syringe	Pirlimycin Hydrochloride (50MG)	Lactating	9 days-meat/ 36 hrs milk
Zoetis	034114	Pirsue Sterile Solution Syringe	144 Count	Box	Pirlimycin Hydrochloride (50MG)	Lactating	9 days-meat/ 36 hrs milk
Zoetis	034115	Spectramast LC	12 Count	Syringe	Ceftiofur Equivalents (125MG), Microcrystalline Wax (700MG), Labrafil M 1944CS (500MG)	Lactating	2 days meat/ 72 hrs milk
Zoetis	034116	Spectramast LC Syringe	144 Count	Box	Ceftiofur Equivalents (125MG), Microcrystalline Wax (700MG), Labrafil M 1944CS (500MG)	Lactating	2 days meat/ 72 hrs milk
Boehringer Ingelheim	034456	Tomorrow Infusion 10mL Tube	144 Count	Pail	Cephapirin Benzathine (300MG)	Dry	42 days-meat/ 72 hrs post calving
Boehringer Ingelheim	034081	Tomorrow Mastitis Treatment 10mL Syringe	12 Count	Box	Cephapirin Benzathine (300MG)	Dry	42 days-meat/ 72 hrs post calving
Boehringer Ingelheim	004122	Dry Clox Infusion	144 SYR	Pail	Cloxacillin Benzathine (500MG)	Dry	30 days-meat/ 30 days pre calving
Merck	021923	Orbenin-DC Dry Cow Intramammary Infusion	12 x 7.5gm	Box	Benazthine Cloxacillin (500MG)	Dry	28 days-meat/ 28 days pre calving
WG Critical Care	050474	Quartermaster Suspension	144 x 10mL	Pail	Procaine Penicillin G (1,000,000 units), Dihydrostreptomycin (1GM)	Dry	60 days- meat/ 96 hrs milk
WG Critical Care	050475	Quartermaster Suspension	12 x 10mL	Box	Procaine Penicillin G (1,000,000 units), Dihydrostreptomycin (1GM)	Dry	60 days- meat/ 96 hrs milk
Zoetis	019993	Albadry Plus Mastitis Tubes	144 x 10mL	Box	Penicillin G Procaine (200,000 IU) Novibiocin Sodium (400MG), Chlorobutanol Anhydrous (50MG)	Dry	30 days-meat/ 72 hrs post calving
Zoetis	034200	Orbeseal Mastitis Tubes	144 Count	Pail	65% Bismuth Subnitrate (4G)	Dry	0 days- meat/ 0 days milk
Zoetis	034117	Spectramast DC	12 Count	Syringe	Ceftiofur Equivalents (500MG), Microcrystalline Wax (700MG), Labrafil M 1944CS (500MG)	Dry	16 days-meat/ 0 days milk
Zoetis	034118	Spectramast DC Syringe	144 Count	Box	Ceftiofur Equivalents (500MG), Microcrystalline Wax (700MG), Labrafil M 1944CS (500MG)	Dry	16 days-meat/ 0 days milk

